

SERVING SCIENTISTS SINCE 1984

DESIGNED AND ASSEMBLED

MICRONIC WARRANTY
OF 1 YEAR

VERSATILE SAMPLE MANAGEMENT EQUIPMENT

Micronic Tube Handler HT700

Next level sorting of your valuable liquid and solid samples









Ordering information

Catalog numberDescriptionUnitMP35500Micronic Tube Handler HT700 (pick and place + 2D reading)1 unit

Superior features

Capacity

- Minimal footprint increase with doubling of capacity
- Source rack capacity: 24
- Target rack capacity: 1-3

Throughput

- Faster operation of 800 tubes per hour due to minimum distance between the sample to be picked and the target rack
- Picking and 2D code reading a high number of samples with a space efficient design
- Features
- 2D code reading
- 1D code reading (source and target rack)
- Easy to navigate and intuitive interface
- Automatic linking between rack ID and rack position
- Loading CSV-file pick list from USB thumb drive
- Extra grouping field in pick list for sample clustering in target rack
- Output CSV-file of rack content
- Enclosed system for operator safety
 Separate safety stop ensures sample/product protection
- Compatible labware
- 96-well format Micronic racks (ANSI/SLAS standard)
- Capped 0.30ml and 0.40ml tubes
 Capped and uncapped 0.50ml, 0.75ml, 0.80ml, 1.10ml and
- 1.40ml Micronic tubes
 Contact us for compatibility with other labware brands
- Operation conditions

 Room temperature operation



