


30 years of experience  
in traceable sample  
storage

Class 7 clean room  
injection molding  
and assembly

Production in  
the Netherlands and  
the United States

## MICRONIC TUBE STORAGE RACK RANGE

### Loborack-96 Holds 96 Individual Tubes



### Superior features

- Based on the Society for Laboratory Automations and Screening (ANSI/SLAS) standards for storage racks
- Standard available with laser etched 1D rack barcode
- Secure and safe locking system
- Twist-lock
- Tapered walls allow precise tube placement into racks (sure shot design)
- Open rack bottom
- Compatible with the Micronic Traxer Code Readers
- Class 7 cleanroom production
- RNase/DNase and Pyrogen free

### Ordering information

Catalog number	Description	Unit
MPW51013BC3	Loborack-96 white, laser etched (including low cover) (barcoded A1-H1 side)	10 racks
MPW51004BC3	Loborack-96 white, laser etched (including standard cover) (barcoded A1-H1 side)	10 racks
MPW51015BC3	Loborack-96 white, laser etched (including high cover) (barcoded A1-H1 side)	10 racks



To contact Micronic, email to [sales@micronic.com](mailto:sales@micronic.com) or call us at:

**Micronic Europe** (Europe, Middle-East and Africa): +31 (0)320 277 070 | **Micronic America** (United States, Latin America and Canada): +1 484 480 3372 | **Micronic Asia** (Asia and Pacific): +31 (0)320 277 099