





MEDIA BOTTLES PRODUCT DATA SHEET												
Bottle Material	PETG, PET											
Cap Material	HDPE											
Temperature Range	PETG -70 to 60°C						PET -70 to 60°C					
Autoclavable	No						No					
Leak proof Tested	715 mmHg pressure at room temperature						715 mmHg pressure at room temperature					
Hazards Data	Health 0 Non-hazardous Flammability 1 Flammable (Class A fire) Reactivity 0 Non-reactive <small>*Consult the polymer MSDS for further information on potential hazards associated with this material.</small>						Health 0 Non-hazardous Flammability 1 Flammable (Class A fire) Reactivity 0 Non-reactive <small>*Consult the polymer MSDS for further information on potential hazards associated with this material.</small>					
Type	PETG							PET				
Catalogue Number	5600030C	5600060C	5600125C	5600250C	5600500C	5601000C	5602000C	5610125C	5610250C	5610500C	5611000C	
Description	CAPP Media Bottle w/ Cap, PETG, 30mL, sterile	CAPP Media Bottle w/ Cap, PETG, 60mL, sterile	CAPP Media Bottle w/ Cap, PETG, 125mL, sterile	CAPP Media Bottle w/ Cap, PETG, 250mL, sterile	CAPP Media Bottle w/ Cap, PETG, 500mL, sterile	CAPP Media Bottle w/ Cap, PETG, 1L, sterile	CAPP Media Bottle w/ Cap, PETG, 2 L, sterile	CAPP Media Bottle w/ Cap, PET, 125mL, sterile	CAPP Media Bottle w/ Cap, PET, 250mL, sterile	CAPP Media Bottle w/ Cap, PET, 500mL, sterile	CAPP Media Bottle w/ Cap, PET, 1L, sterile	
Sales Unit	Case / 4 x 24 pcs.	Case / 4 x 24 pcs.	Case / 2 x 24 pcs.	Case / 2 x 24 pcs.	Case / 2 x 12 pcs.	Case / 2 x 12 pcs.	Case / 2 x 6 pcs.	Case / 4 x 24 pcs.	Case / 2 x 30 pcs.	Case / 2 x 24 pcs.	Case / 2 x 12 pcs.	
GTIN	4255662905019	4255662905026	4255662905033	4255662905040	4255662905057	4255662905064	4255662905071	4255662905088	4255662905095	4255662905101	4255662905118	
Sterile	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Volume (ml)	30	60	125	250	500	1000	2000	125	250	500	1000	
Image												
B (With closure) mm	63,3	82,6	107,6	144,6	179,0	220,6	271,5	107,4	146,4	176,7	217,9	
C (Without Closure) mm	61,6	81,2	103,7	140,7	175,0	216,7	265,5	102,7	142,5	172,8	214,0	
D Square (mm)	39,0	41,2	53,4	61,0	74,6	93,1	116,7	52,3	60,4	77,6	92,2	
A Bottle neck (ID)	14,0	18,0	28,5	28,5	28,5	28,5	39,75	28,5	28,5	28,5	28,5	
Closure Size (mm)	20-415	24-415	38-430	38-430	38-430	38-430	53B	38-430	38-430	38-430	38-430	
Total weight (g)	13,3	20,0	44,0	56,0	111,0	182,0	337	39,1	45,0	72,0	102,5	
Production Standards and Procedures	Product is produced under controlled environment. In the whole assembly area caps, masks, gloves and gowns are worn. Lot size is variable due to continuous molding process.											



MEDIA BOTTLES											
PRODUCT DATA SHEET											
Bottle Material	PETG, PET										
Cap Material	HDPE										
Temperature Range	PETG -70 to 60°C						PET -70 to 60°C				
Autoclavable	No						No				
Leak proof Tested	715 mmHg pressure at room temperature						715 mmHg pressure at room temperature				
Hazards Data	Health 0 Non-hazardous Flammability 1 Flammable (Class A fire) Reactivity 0 Non-reactive <small>*Consult the polymer MSDS for further information on potential hazards associated with this material.</small>						Health 0 Non-hazardous Flammability 1 Flammable (Class A fire) Reactivity 0 Non-reactive <small>*Consult the polymer MSDS for further information on potential hazards associated with this material.</small>				
Type	PETG						PET				
Catalogue Number	5500030C	5500060C	5500125C	5500250C	5500500C	5501000C	5502000C	5510125C	5510250C	5510500C	5511000C
Description	CAPP Media Bottle w/ Cap, PETG, 30mL, Bulk, non-sterile	CAPP Media Bottle w/ Cap, PETG, 60mL, Bulk, non-sterile	CAPP Media Bottle w/ Cap, PETG, 125mL, Bulk, non-sterile	CAPP Media Bottle w/ Cap, PETG, 250mL, Bulk, non-sterile	CAPP Media Bottle w/ Cap, PETG, 500mL, Bulk, non-sterile	CAPP Media Bottle w/ Cap, PETG, 1L, Bulk, non-sterile	CAPP Media Bottle w/ Cap, PETG, 2L, Bulk, non-sterile	CAPP Media Bottle w/ Cap, PET, 125mL, non-sterile	CAPP Media Bottle w/ Cap, PET, 250mL, non-sterile	CAPP Media Bottle w/ Cap, PET, 500mL, non-sterile	CAPP Media Bottle w/ Cap, PET, 1L, non-sterile
Sales Unit	Case / 12 x 24 pcs.	Case / 8 x 24 pcs.	Case / 4 x 24 pcs.	Case / 2 x 30 pcs.	Case / 2 x 24 pcs.	Case / 2 x 12 pcs.	Case / 2 x 6 pcs.	Case / 4 x 24 pcs.	Case / 2 x 30 pcs.	Case / 2 x 24 pcs.	Case / 2 x 12 pcs.
GTIN	4255662905682	4255662905699	4255662905705	4255662905712	4255662905729	4255662905736	4255662905743	4255662905750	4255662905767	4255662905774	4255662905781
Sterile											
Volume (ml)	30	60	125	250	500	1000	2000	125	250	500	1000
Image											
B (With closure) mm	63,3	82,6	107,6	144,6	179,0	220,6	271,5	107,4	146,4	176,7	217,9
C (Without Closure) mm	61,6	81,2	103,7	140,7	175,0	216,7	265,5	102,7	142,5	172,8	214,0
D Square (mm)	39,0	41,2	53,4	61,0	74,6	93,1	116,7	52,3	60,4	77,6	92,2
A Bottle neck (ID)	14,0	18,0	28,5	28,5	28,5	28,5	39,75	28,5	28,5	28,5	28,5
Closure Size (mm)	20-415	24-415	38-430	38-430	38-430	38-430	53B	38-430	38-430	38-430	38-430
Total weight (g)	13,3	20,0	44,0	56,0	111,0	182,0	337	39,1	45,0	72,0	102,5
Production Standards and Procedures	Product is produced under controlled environment. In the whole assembly area caps, masks, gloves and gowns are worn. Lot size is variable due to continuous molding process.										



MEDIA BOTTLES PRODUCT DATA SHEET					
Bottle Material	PETG, PET				
Cap Material	HDPE				
Temperature Range	PETG -70 to 60°C				
Autoclavable	No				
Leak proof Tested	715 mmHg pressure at room temprature				
Hazards Data	Health 0 Non-hazardous Flammability 1 Flammable (Class A fire) Reactivity 0 Non-reactive <small>*Consult the polymer MSDS for further information on potential hazards associated with this material.</small>				
Type	PETG				
Catalogue Number	5620030C	5620060C	5620125C	5620250C	5620500C
Description	CAPP Media Bottle w/ Cap, PETG, 30mL, Bulk, sterile	CAPP Media Bottle w/ Cap, PETG, 60mL, Bulk, sterile	CAPP Media Bottle w/ Cap, PETG, 125mL, Bulk, sterile	CAPP Media Bottle w/ Cap, PETG, 250mL, Bulk, sterile	CAPP Media Bottle w/ Cap, PETG, 500mL, Bulk, sterile
Sales Unit	Case / 12 x 24 pcs.	Case / 8 x 24 pcs.	Case / 4 x 24 pcs.	Case / 2 x 30 pcs.	Case / 2 x 24 pcs.
GTIN	4255662905798	4255662905804	4255662905811	4255662905828	4255662905835
Sterile	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Volume (ml)	30	60	125	250	500
Image					
B (With closure) mm	63,3	82,6	107,6	144,6	179,0
C (Without Closure) mm	61,6	81,2	103,7	140,7	175,0
D Square (mm)	39,0	41,2	53,4	61,0	74,6
A Bottle neck (ID)	14,0	18,0	28,5	28,5	28,5
Closure Size (mm)	20-415	24-415	38-430	38-430	38-430
Total weight (g)	13,3	20,0	44,0	56,0	111,0
Production Standards and Procedures	Product is produced under controlled environment. In the whole assembly area caps, masks, gloves and gowns are worn. Lot size is variable due to continuous molding process.				



Technical Drawing