



# TECHNICAL SERVICE PLANS

## Scinomix Sci-Print VXL



|                                                                               | LEVEL ONE | LEVEL TWO | LEVEL THREE |
|-------------------------------------------------------------------------------|-----------|-----------|-------------|
| One scheduled preventative maintenance (PM) inspection per instrument         | ✓         | ✓         | ✓           |
| One on-site service visit for emergency needs                                 |           | ✓         | ✓           |
| One additional on-site service visit for emergency needs or second PM         |           |           | ✓           |
| Travel expenses, including overnight stays                                    | ✓         | ✓         | ✓           |
| Performance tests of included equipment and parts                             | ✓         | ✓         | ✓           |
| Thorough cleaning and lubrications of rails and axis during the PM visit      | ✓         | ✓         | ✓           |
| Functionality testing of sensors, motors, and accessories during the PM visit | ✓         | ✓         | ✓           |
| Inspection and when required the replacement of parts*                        | ✓         | ✓         | ✓           |
| Scinomix control software upgrades                                            |           |           | ✓           |

\*list of included parts is available on the second page of this flyer

**HAVE ANY QUESTIONS  
ABOUT OUR SERVICES?  
REACH OUT:**

**NBS Scientific FR**  
**info@nbsscientific.fr**  
**+33 (0)478 964 700**  
**www.nbsscientific.fr**

**OUR TEAM IS  
HERE TO HELP!**

## **TECHNICAL SERVICE PLANS** **Scinomix Sci-Print VXL**

### **INCLUDED SPARE PARTS:**

| <b>Part Number</b> | <b>Part Description</b>                                                                              |
|--------------------|------------------------------------------------------------------------------------------------------|
| SC2000-090         | Drive Print Roller for A2+ and SQUIX                                                                 |
| SC2000-019         | Bearing, 1/4" Bore x 3/8" OD, Flanged, Shielded with Extended Inner Ring, R168-2Z                    |
| SC2000-274         | M3-0.5 x 10mm L, SS FHCS                                                                             |
| SC2000-091         | Rewind Assist Roller for A2+ and SQUIX                                                               |
| SC2000-072         | Bearing, 3/16" Bore x 1/2" OD, Shielded with Extended Inner Ring, R3-27                              |
| SC2000-284         | PM Stickers                                                                                          |
| SC2000-1410        | Bearing, 5/16" Bore x 1/2" OD, Flanged, Shielded with Extended Inner Ring, R1810-2Z                  |
| SC2000-0913        | Bearing, Ball, Flanged Shielded, Extended Inner Ring, 1/4 ID x 1/2 OD x 3/16 W                       |
| SC2000-2207        | 18-8 Stainless Steel Shoulder Screw 1/4" Shoulder Diameter, 7/16" Shoulder Length, 10-24 Thread Size |
| SC2000-019         | Bearing, 1/4" Bore x 3/8" OD, Flanged, Shielded with Extended Inner Ring, R168-2Z                    |
| SC2000-628         | 1/8 x 9in Spindle Drive Round Belt                                                                   |
| SC2000-1053        | Retaining Ring: External: 0.250in Shaft                                                              |

**NBS Scientific FR**  
**93 Avenue du Chateau**  
**13940, Mollégès**  
**France**  
**info@nbsscientific.fr**  
**+33 (0)478 964 700**

**TECHNICAL SERVICE**  
**Meerpaalweg 6**  
**1332 BB, Almere**  
**The Netherlands**  
**+31 (0)36 549 1010**  
**ts@nbsscientific.nl**